


The Autoprep



Stainless Steel Case

method of preparing peripheral blood films using the wedge technique.

Autoprep features include:

User Friendly Design • Simple Setup & Use • Easy Clean Up • Low Maintenance • Easy to Follow Operator Manual • Constant Dispense Rate • Pause Control • Smear Control Knob • Standard Wedge Configuration • Fresh & Anti-Coagulated Blood • Stationary or Portable use • Stainless Steel Case • Durable Construction • 1 year warranty

Specifications: 6' x 4' x 6' (L x W x H) 3.2 lbs. • 15cm x 10cm x 15cm (L x W x H) (1.45 Kg)




Step 1: Place slide



Step 2: Drop blood



Step 3 : Press down



Step 4: Release



Step 5 : Clean glass

www.SedonaLab.com
Info@Sedonalab.com